

# Constant - Temperature Drying Oven

1. External material: excellent cold-rolled steel with static coating.
2. Stainless steel inner chamber;
3. PID micro-computer intelligent temperature controller, over-temperature protection; LED display;
4. Double-layer glass observation window, adjustable airtight buckle lock.
5. No blowing, natural convection heating circulation.

Optional:

1. LCD digital screen
2. Independent temperature controller



Model	High brightness digital tube	BOV-T30CII	BOV-T45C	BOV-V65C	BOV-V125C
voltage		220V/50Hz			
Tem. Range(°C)		RT+5~300°C			
Tem. Motion(°C)		±1°C			
Tem.Resolution Ratio(°C)		0.1°C			
Ambient tem.		5-40°C			
Power(Kw)		0.8	1	1.5	2.5
Inner Champer Size(mm)		310×310×310	350×350×350	400×360×450	500×450×550
Exterior Size(mm)		450×490×685	490×530×730	540×540×830	640×630×930
Packing Size (mm)		545×520×750	590×560×790	610×610×890	710×700×990
Load per rack(Kg)		15/kg 2 for standard configuration			
N.W.(Kg)		33	37	44	60
G.W. (Kg)		37	43	49	66
Timing		1-9999min			